




PPP2R1B Antibody

Product Code	CSB-PA003833
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30154
Immunogen	Synthesized peptide derived from the C-terminal region of Human PP2A-A β .
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IHC:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP2R1B; Serine/threonine-protein phosphatase 2A 65 kDa regulatory subunit A beta isoform; PP2A subunit A isoform PR65-beta; PP2A subunit A isoform R1-beta
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2R1B

Image



Western Blot analysis of Jurkat cells using PPP2A-A β Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.