





PPP2CA Antibody

Product Code	CSB-PA003834
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P67775
Immunogen	Synthesized peptide derived from Human PP2A-C α around the non-phosphorylation site of Y307.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP2CA; Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform; PP2A-alpha; Replication protein C; RP-C
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2CA

Image



Western Blot analysis of A549 cells using PPP2A-C α Polyclonal Antibody



Western Blot analysis of Jurkat cells using PP2A-C α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.