




# PPM1K Antibody

|                            |                                                                                                                                                                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA003835                                                                                                                                                                                                                                    |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                   |
| <b>Uniprot No.</b>         | Q8N3J5                                                                                                                                                                                                                                          |
| <b>Immunogen</b>           | Synthesized peptide derived from the Internal region of Human PP2Cκ.                                                                                                                                                                            |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                          |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                           |
| <b>Tested Applications</b> | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000                                                                                                                                                                                 |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                          |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                                                                                                                  |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                         |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                           |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                             |
| <b>Alias</b>               | PPM1K; PP2CM; Protein phosphatase 1K; mitochondrial; PP2C domain-containing protein phosphatase 1K; PP2C-like mitochondrial protein; PP2C-type mitochondrial phosphoprotein phosphatase; PTMP; Protein phosphatase 2C isoform kappa; PP2C-kappa |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                            |
| <b>Target Names</b>        | PPM1K                                                                                                                                                                                                                                           |

## Image



Western Blot analysis of 293 cells using PP2Cκ Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.