




PPP1R14D Antibody

Product Code	CSB-PA003840
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXH3
Immunogen	Synthesized peptide derived from the Internal region of Human PPP1R14D.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PPP1R14D; GBPI; Protein phosphatase 1 regulatory subunit 14D; Gastrointestinal and brain-specific PP1-inhibitory protein 1; GBPI-1
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R14D

Image



Western Blot analysis of COLO205 cells using PPP1R14D Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.