




PGR Antibody

Product Code	CSB-PA003846
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06401
Immunogen	Synthesized peptide derived from Human PR around the non-phosphorylation site of S400.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
Immunogen Species	Homo sapiens (Human)
Target Names	PGR

Image



Western Blot analysis of HeLa cells using PR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.