




PTPN1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003883 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18031 |
| Immunogen | Synthesized peptide derived from Human PTP1B around the non-phosphorylation site of S50. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat,Monkey |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | PTPN1; PTP1B; Tyrosine-protein phosphatase non-receptor type 1; Protein-tyrosine phosphatase 1B; PTP-1B |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PTPN1 |

Image



Western Blot analysis of COS7 cells using PTP1B Polyclonal Antibody