




PTTG1 Antibody

Product Code	CSB-PA003885
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95997
Immunogen	Synthesized peptide derived from the Internal region of Human PTTG1/2/3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PTTG1; EAP1; PTTG; TUTR1; Securin; Esp1-associated protein; Pituitary tumor-transforming gene 1 protein; Tumor-transforming protein 1; hPTTG
Immunogen Species	Homo sapiens (Human)
Target Names	PTTG1

Image



Western Blot analysis of 293 cells using PTTG1/2/3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.