




BBC3 Antibody

Product Code	CSB-PA003890
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXH1
Immunogen	Synthesized peptide derived from the C-terminal region of Human PUMA.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BBC3; PUMA; Bcl-2-binding component 3; JFY-1; p53 up-regulated modulator of apoptosis
Immunogen Species	Homo sapiens (Human)
Target Names	BBC3

Image



Western Blot analysis of 293 cells using PUMA Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.