





RAB34 Antibody

Product Code	CSB-PA003895
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZG1
Immunogen	Synthesized peptide derived from the C-terminal region of Human Rab 34.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAB34; RAB39; RAH; Ras-related protein Rab-34; Ras-related protein Rab-39; Ras-related protein Rah
Immunogen Species	Homo sapiens (Human)
Target Names	RAB34

Image



Western Blot analysis of 293 cells using Rab 34 Polyclonal Antibody



Western blot analysis of 293T lysis using Rab 34 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.