




RRAD Antibody

Product Code	CSB-PA003910
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55042
Immunogen	Synthesized peptide derived from the Internal region of Human Rad GTPase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RRAD; RAD; GTP-binding protein RAD; RAD1; Ras associated with diabetes
Immunogen Species	Homo sapiens (Human)
Target Names	RRAD

Image



Western Blot analysis of MCF7 cells using Rad GTPase Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.