




RAD17 Antibody

Product Code	CSB-PA003912
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75943
Immunogen	Synthesized peptide derived from Human Rad17 around the non-phosphorylation site of S645.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAD17; R24L; Cell cycle checkpoint protein RAD17; hRad17; RF-C/activator 1 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	RAD17

Image



Western Blot analysis of HuvEc cells using Rad17 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.