





RAD23B Antibody

Product Code	CSB-PA003914
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54727
Immunogen	Synthesized peptide derived from the N-terminal region of Human Rad23B.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAD23B; UV excision repair protein RAD23 homolog B; HR23B; hHR23B; XP-C repair-complementing complex 58 kDa protein; p58
Immunogen Species	Homo sapiens (Human)
Target Names	RAD23B

Image



Western Blot analysis of HuvEc cells using Rad23B Polyclonal Antibody



Western Blot analysis of Jurkat K562 cells using Rad23B Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.