





RAD52 Antibody

Product Code	CSB-PA003918
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43351
Immunogen	Synthesized peptide derived from Human Rad52 around the non-phosphorylation site of Y104.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAD52; DNA repair protein RAD52 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	RAD52

Image



Western Blot analysis of 293 cells using Rad52 Polyclonal Antibody



Western Blot analysis of Jurkat cells using Rad52 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.