




RRAGC Antibody

Product Code	CSB-PA003932
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB90
Immunogen	Synthesized peptide derived from the C-terminal region of Human Rag C.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RRAGC; Ras-related GTP-binding protein C; Rag C; RagC; GTPase-interacting protein 2; TIB929
Immunogen Species	Homo sapiens (Human)
Target Names	RRAGC

Image



Western Blot analysis of LOVO cells using Rag C Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.