




RANBP6 Antibody

Product Code	CSB-PA003936
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60518
Immunogen	Synthesized peptide derived from the Internal region of Human Ran BP-6.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RANBP6; Ran-binding protein 6; RanBP6
Immunogen Species	Homo sapiens (Human)
Target Names	RANBP6

Image



Western Blot analysis of COLO205 cells using Ran BP-6 Polyclonal Antibody