





RAP1GAP Antibody

Product Code	CSB-PA003941
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47736
Immunogen	Synthesized peptide derived from the C-terminal region of Human Rap1GAP.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RAP1GAP; KIAA0474; RAP1GA1; Rap1 GTPase-activating protein 1; Rap1GAP; Rap1GAP1
Immunogen Species	Homo sapiens (Human)
Target Names	RAP1GAP

Image



Western Blot analysis of COLO205 cells using Rap1GAP Polyclonal Antibody



Western Blot analysis of PC3 NIH-3T3 cells using Rap1GAP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.