




RASSF6 Antibody

Product Code	CSB-PA003945
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZTQ3
Immunogen	Synthesized peptide derived from the Internal region of Human RASSF6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RASSF6; Ras association domain-containing protein 6
Immunogen Species	Homo sapiens (Human)
Target Names	RASSF6

Image



Western Blot analysis of HeLa cells using RASSF6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.