




RBAK Antibody

Product Code	CSB-PA003952
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYW8
Immunogen	Synthesized peptide derived from the N-terminal region of Human RBAK.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RBAK; ZNF769; RB-associated KRAB zinc finger protein; RB-associated KRAB repressor; hRBaK; Zinc finger protein 769
Immunogen Species	Homo sapiens (Human)
Target Names	RBAK

Image



Western Blot analysis of HepG2 cells using RBAK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.