




RFX3 Antibody

Product Code	CSB-PA003966
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48380
Immunogen	Synthesized peptide derived from the C-terminal region of Human RFX3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RFX3; Transcription factor RFX3; Regulatory factor X 3
Immunogen Species	Homo sapiens (Human)
Target Names	RFX3

Image



Western Blot analysis of RAW264.7 cells using RFX3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.