




RGS10 Antibody

Product Code	CSB-PA003969
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43665
Immunogen	Synthesized peptide derived from the Internal region of Human RGS10.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RGS10; Regulator of G-protein signaling 10; RGS10
Immunogen Species	Homo sapiens (Human)
Target Names	RGS10

Image



Western Blot analysis of Jurkat cells using RGS10 Polyclonal Antibody