




RPL10L Antibody

Product Code	CSB-PA003978
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96L21
Immunogen	Synthesized peptide derived from the C-terminal region of Human Ribosomal Protein L10L.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPL10L; 60S ribosomal protein L10-like
Immunogen Species	Homo sapiens (Human)
Target Names	RPL10L

Image



Western Blot analysis of HeLa cells using Ribosomal Protein L10L Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.