





RPS4X Antibody

Product Code	CSB-PA004003
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62701
Immunogen	Synthesized peptide derived from the Internal region of Human Ribosomal Protein S4X.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS4X; CCG2; RPS4; SCAR; 40S ribosomal protein S4; X isoform; SCR10; Single copy abundant mRNA protein
Immunogen Species	Homo sapiens (Human)
Target Names	RPS4X

Image



Western Blot analysis of COLO205 cells using Ribosomal Protein S4X Polyclonal Antibody



Western Blot analysis of customer's (cat sample) using Ribosomal Protein S4X Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.