





RIMS4 Antibody

Product Code	CSB-PA004009
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H426
Immunogen	Synthesized peptide derived from the Internal region of Human Rim4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RIMS4; C20orf190; Regulating synaptic membrane exocytosis protein 4; RIM4 gamma; Rab3-interacting molecule 4; RIM 4
Immunogen Species	Homo sapiens (Human)
Target Names	RIMS4

Image



Western Blot analysis of HeLa cells using Rim4 Polyclonal Antibody



Western blot analysis of 3T3 lysis using Rim4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.