





MPRIP Antibody

Product Code	CSB-PA004012
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6WCQ1
Immunogen	Synthesized peptide derived from the Internal region of Human M-RIP.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MPRIP; KIAA0864; MRIP; RHOIP3; Myosin phosphatase Rho-interacting protein; M-RIP; Rho-interacting protein 3; RIP3; p116Rip
Immunogen Species	Homo sapiens (Human)
Target Names	MPRIP

Image



Western Blot analysis of 293 cells using M-RIP Polyclonal Antibody



Western Blot analysis of K562 Jurkat RAT-MUSCLE cells using M-RIP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.