




ROCK2 Antibody

Product Code	CSB-PA004018
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75116
Immunogen	Synthesized peptide derived from Human Rock-2 around the non-phosphorylation site of Y722.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ROCK2; KIAA0619; Rho-associated protein kinase 2; Rho kinase 2; Rho-associated; coiled-coil-containing protein kinase 2; Rho-associated, coiled-coil-containing protein kinase II; ROCK-II; p164 ROCK-2
Immunogen Species	Homo sapiens (Human)
Target Names	ROCK2

Image



Western Blot analysis of RAT-MUSCLE cells using Rock-2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.