




RPS6KA6 Antibody

Product Code	CSB-PA004037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK32
Immunogen	Synthesized peptide derived from the C-terminal region of Human Rsk-4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RPS6KA6; RSK4; Ribosomal protein S6 kinase alpha-6; S6K-alpha-6; 90 kDa ribosomal protein S6 kinase 6; p90-RSK 6; p90RSK6; Ribosomal S6 kinase 4; RSK-4; pp90RSK4
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KA6

Image



Western Blot analysis of 453 293T AD293 HELA HepG2 cells using Rsk-4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.