




SHC1 Antibody

Product Code	CSB-PA004082
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29353
Immunogen	Synthesized peptide derived from Human Shc around the non-phosphorylation site of Y427.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SHC1; SHC; SHCA; SHC-transforming protein 1; SHC-transforming protein 3; SHC-transforming protein A; Src homology 2 domain-containing-transforming protein C1; SH2 domain protein C1
Immunogen Species	Homo sapiens (Human)
Target Names	SHC1

Image



Western Blot analysis of HepG2 cells using Shc Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.