




SMAD1/SMAD5/SMAD9 Antibody

Product Code	CSB-PA004105
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15797/Q99717/O15198
Immunogen	Synthesized peptide derived from the N-terminal region of Human Smad1/5/9.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SMAD1; BSP1; MADH1; MADR1; Mothers against decapentaplegic homolog 1; MAD homolog 1; Mothers against DPP homolog 1; JV4-1; Mad-related protein 1; SMAD family member 1; SMAD 1; Smad1; hSMAD1; Transforming growth factor-beta-signaling protein
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD1/SMAD5/SMAD9

Image



Western Blot analysis of HeLa cells using Smad1/5/9 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.