





PLK2 Antibody

Product Code	CSB-PA004124
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYY3
Immunogen	Synthesized peptide derived from the Internal region of Human Snk.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLK2; SNK; Serine/threonine-protein kinase PLK2; Polo-like kinase 2; PLK-2; hPlk2; Serine/threonine-protein kinase SNK; hSNK; Serum-inducible kinase
Immunogen Species	Homo sapiens (Human)
Target Names	PLK2

Image



Western Blot analysis of MCF-7 cells using Snk Polyclonal Antibody



Western Blot analysis of HeLa cells using Snk Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.