




EFTUD2 Antibody

Product Code	CSB-PA004127
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15029
Immunogen	Synthesized peptide derived from the Internal region of Human Snrp116.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EFTUD2; KIAA0031; SNRP116; 116 kDa U5 small nuclear ribonucleoprotein component; Elongation factor Tu GTP-binding domain-containing protein 2; SNU114 homolog; hSNU114; U5 snRNP-specific protein; 116 kDa; U5-116 kDa
Immunogen Species	Homo sapiens (Human)
Target Names	EFTUD2

Image



Western Blot analysis of HuvEc cells using Snrp116 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.