





SST Antibody

Product Code	CSB-PA004132
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61278
Immunogen	Synthesized peptide derived from the N-terminal region of Human Somatostatin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SST; Somatostatin; Growth hormone release-inhibiting factor
Immunogen Species	Homo sapiens (Human)
Target Names	SST

Image



Western Blot analysis of A549 cells using Somatostatin Polyclonal Antibody



Western Blot analysis of PC12 cells using Somatostatin Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.