




SUPT3H Antibody

Product Code	CSB-PA004141
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75486
Immunogen	Synthesized peptide derived from the Internal region of Human SPT3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SUPT3H; SPT3; Transcription initiation protein SPT3 homolog; SPT3-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	SUPT3H

Image



Western Blot analysis of HepG2 cells using SPT3 Polyclonal Antibody