




SUPT3H Antibody

Product Code	CSB-PA004141
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75486
Immunogen	Synthesized peptide derived from the Internal region of Human SPT3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SUPT3H; SPT3; Transcription initiation protein SPT3 homolog; SPT3-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	SUPT3H

Image



Western Blot analysis of HepG2 cells using SPT3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.