





SREBF1 Antibody

Product Code	CSB-PA004153
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36956
Immunogen	Synthesized peptide derived from Human SREBP-1 around the non-phosphorylation site of S439.
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SREBF1; BHLHD1; SREBP1; Sterol regulatory element-binding protein 1; SREBP-1; Class D basic helix-loop-helix protein 1; bHLHd1; Sterol regulatory element-binding transcription factor 1
Immunogen Species	Homo sapiens (Human)
Target Names	SREBF1

Image



Western Blot analysis of hela cells using SREBP-1 Polyclonal Antibody



Western blot analysis of 3T3 lysis using SREBP-1 antibody.