




ZSCAN26 Antibody

Product Code	CSB-PA004154
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16670
Immunogen	Synthesized peptide derived from the Internal region of Human SRE-ZBP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF187; ZSCAN26; Zinc finger protein 187; Protein SRE-ZBP; Zinc finger and SCAN domain-containing protein 26
Immunogen Species	Homo sapiens (Human)
Target Names	ZSCAN26

Image



Western Blot analysis of Jurkat cells using SRE-ZBP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.