




SRSF4 Antibody

Product Code	CSB-PA004160
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08170
Immunogen	Synthesized peptide derived from the Internal region of Human SRp75.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SRSF4; SFRS4; SRP75; Serine/arginine-rich splicing factor 4; Pre-mRNA-splicing factor SRP75; SRP001LB; Splicing factor; arginine/serine-rich 4
Immunogen Species	Homo sapiens (Human)
Target Names	SRSF4

Image



Western Blot analysis of NIH-3T3 cells using SRp75 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.