




ST18 Antibody

Product Code	CSB-PA004164
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60284
Immunogen	Synthesized peptide derived from the N-terminal region of Human ST18.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ST18; KIAA0535; ZNF387; Suppression of tumorigenicity 18 protein; Zinc finger protein 387
Immunogen Species	Homo sapiens (Human)
Target Names	ST18

Image



Western Blot analysis of HuvEc cells using ST18 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.