





STAT5A Antibody

Product Code	CSB-PA004179
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42229
Immunogen	Synthesized peptide derived from the Internal region of Human Stat5a.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STAT5A; STAT5; Signal transducer and activator of transcription 5A
Immunogen Species	Homo sapiens (Human)
Target Names	STAT5A

Image



Western Blot analysis of 293 cells using Stat5a Polyclonal Antibody



Western Blot analysis of K562 cells using Stat5a Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.