





STK33 Antibody

Product Code	CSB-PA004188
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYT3
Immunogen	Synthesized peptide derived from the N-terminal region of Human STK33.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	STK33; Serine/threonine-protein kinase 33
Immunogen Species	Homo sapiens (Human)
Target Names	STK33

Image



Western Blot analysis of 293 cells using STK33 Polyclonal Antibody



Western Blot analysis of HuvEc Jurkat cells using STK33 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.