




SYT13 Antibody

Product Code	CSB-PA004206
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8C5
Immunogen	Synthesized peptide derived from the Internal region of Human Synaptotagmin XIII.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SYT13; KIAA1427; Synaptotagmin-13; Synaptotagmin XIII; SytXIII
Immunogen Species	Homo sapiens (Human)
Target Names	SYT13

Image



Western Blot analysis of HuvEc cells using Synaptotagmin XIII Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.