




SNCAIP Antibody

Product Code	CSB-PA004207
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6H5
Immunogen	Synthesized peptide derived from the C-terminal region of Human Synphilin-1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SNCAIP; Synphilin-1; Sph1; Alpha-synuclein-interacting protein
Immunogen Species	Homo sapiens (Human)
Target Names	SNCAIP

Image



Western Blot analysis of COLO205 cells using Synphilin-1 Polyclonal Antibody