




TAS2R39 Antibody

Product Code	CSB-PA004217
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59534
Immunogen	Synthesized peptide derived from the Internal region of Human T2R39.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TAS2R39; Taste receptor type 2 member 39; T2R39; Taste receptor type 2 member 57; T2R57
Immunogen Species	Homo sapiens (Human)
Target Names	TAS2R39

Image



Western Blot analysis of HepG2 cells using TAS2R39 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.