




TNKS1BP1 Antibody

Product Code	CSB-PA004221
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C2
Immunogen	Synthesized peptide derived from the C-terminal region of Human TAB182.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNKS1BP1; KIAA1741; TAB182; 182 kDa tankyrase-1-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	TNKS1BP1

Image



Western Blot analysis of 293 cells using TAB182 Polyclonal Antibody