




TNKS1BP1 Antibody

Product Code	CSB-PA004221
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C2
Immunogen	Synthesized peptide derived from the C-terminal region of Human TAB182.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNKS1BP1; KIAA1741; TAB182; 182 kDa tankyrase-1-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	TNKS1BP1

Image



Western Blot analysis of 293 cells using TAB182 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.