




TADA2A Antibody

Product Code	CSB-PA004223
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75478
Immunogen	Synthesized peptide derived from the Internal region of Human TADA2L.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TADA2A; TADA2L; KL04P; Transcriptional adapter 2-alpha; Transcriptional adapter 2-like; ADA2-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	TADA2A

Image



Western Blot analysis of COLO205 cells using TADA2L Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.