




TAF5L Antibody

Product Code	CSB-PA004226
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75529
Immunogen	Synthesized peptide derived from the Internal region of Human TAF5L.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TAF5L; PAF65B; TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L; PCAF-associated factor 65 beta; PAF65-beta
Immunogen Species	Homo sapiens (Human)
Target Names	TAF5L

Image



Western Blot analysis of COS7 cells using TAF5L Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.