




MAP3K7 Antibody

Product Code	CSB-PA004227
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43318
Immunogen	Synthesized peptide derived from Human Tak1 around the non-phosphorylation site of T187.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP3K7; TAK1; Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1; TGF-beta-activated kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K7

Image



Western Blot analysis of NIH-3T3 MCF7 cells using Tak1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.