




TBX1 Antibody

Product Code	CSB-PA004238
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43435
Immunogen	Synthesized peptide derived from the C-terminal region of Human TBX1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TBX1; T-box transcription factor TBX1; T-box protein 1; Testis-specific T-box protein
Immunogen Species	Homo sapiens (Human)
Target Names	TBX1

Image



Western Blot analysis of COLO205 cells using TBX1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.