




F3 Antibody

Product Code	CSB-PA004263
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13726
Immunogen	Synthesized peptide derived from Human TF around the non-phosphorylation site of S290.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142
Immunogen Species	Homo sapiens (Human)
Target Names	F3

Image



Western Blot analysis of COS7 cells using TF Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.