





GTF3C3 Antibody

Product Code	CSB-PA004269
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Q9
Immunogen	Synthesized peptide derived from the N-terminal region of Human TFIIC102.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GTF3C3; General transcription factor 3C polypeptide 3; Transcription factor IIIC 102 kDa subunit; TFIIC 102 kDa subunit; TFIIC102; Transcription factor IIIC subunit gamma; TF3C-gamma
Immunogen Species	Homo sapiens (Human)
Target Names	GTF3C3

Image



Western Blot analysis of 293 cells using TFIIC102 Polyclonal Antibody



Western Blot analysis of LOVO NIH-3T3 cells using TFIIC102 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.